**DENTSLEEVE MANOMETRIC CATHETERS  
INSTRUCTIONS FOR USE**

**INTENDED PURPOSE**

Dentsleeve Silicone Manometric Catheters are to be used in conjunction with a manometric pump and a computerized data processing system, to measure muscular contractions along the gastrointestinal system. This manometric catheter is designed to be reusable, and is to be used only under the supervision of a physician who has received adequate training in gastro-intestinal manometry.

**INDICATIONS FOR USE**

Use of this catheter is indicated when measurements of gastro-intestinal tract pressures are judged to be useful for determining management of patients with proven or suspected gastro-intestinal motor disorders.

**CONTRAINDICATIONS**

Where there is a significantly increased risk of pulmonary aspiration associated with intubation which cannot be adequately reduced by practical measures such as a longer than usual period of fasting prior to intubation. The presence of any stenosis, structural deformity or disease that results in an unacceptable increased risk of aspiration, perforation or any other adverse consequence arising from use of this assembly.

**WARNINGS**

This device is supplied non-sterile. Proper aseptic technique and universal precautions (UBP) must apply. Post-procedure manometric catheters are contaminated with body fluid. Immediately clean and disinfect according to the procedures stipulated in the accompanying Instructions for Use.

**INSTRUCTIONS FOR USE**

- Inspect catheter prior to use; do not use if damaged.
  - Water fill channels that will not be used for measurements and seal catheter luer connector with the luer plugs provided.
  - Labeled luer connector ends are to be connected to the corresponding channel on the manometric pump.
  - Post-procedure catheters are contaminated with body fluid; only re-use catheters that have been cleaned and reprocessed according to instructions below.
  - Proper aseptic technique and universal barrier precautions (UBP) must apply.
- POSITIONING OF SLEEVE CATHETERS (IF APPLICABLE)**
- The Sleeve should be positioned so that it remains within the sphincter during its normal movement.
  - After deeper insertion of the Sleeve, pull the catheter back into the sphincter, observing tracings from the Sleeve and side holes at either end of the Sleeve.
  - Recording of a sphincter pattern of pressure from a side hole at either end of the Sleeve indicates placement of the catheter relative to the sphincter.

**CLEANING AND STERILIZATION**

- This catheter is a re-usable device which has been validated as being able to withstand 50 cycles of steam autoclaving.
- If catheter has balloon attached, plug the end of the luer connector to ensure no water enters into the balloon.
  - If catheter has impedance rings, do not immerse the Redel Connector in any liquid.
  - Immediately after extubation, immerse the catheter in a bowl of warm, mild enzymatic detergent and wipe externally several times. Do not allow the uncleaned catheter to dry out. Flush the detergent solution gently through each channel – Do not over-pressurize.
  - Rinse the outside of the catheter in a bowl of clean water. Place the catheter on a clean towel, and cover it with a fold of the towel.
  - Flush each channel first with water, then air – do not over-pressurize. Remove the luer connector plug from the balloon channel and autoclave with the catheter.
  - Package the catheter for autoclaving, and autoclave without delay at 134°C for five minutes at 206kPa or 30psi with a total cycle time of 30 minutes to allow for warm up and cool down time.

**STORAGE**

Ensure all channels and external surfaces are completely dry; store catheter dry in open air. **Do not store in a sealed container.**

Temperature limit:  
15°C - 30°C (60.8°F - 86.6°F)

Humidity limit:  
20% - 75%

**Lifetime of product – 5 years or 50 uses whichever occurs first.**

**Manufacturer's Warranty – 6 months from Manufactured Date for any manufacturer's defects.**

User is to contact the manufacturer or competent authority to report serious incidents involving this device.

Dispose of the device safely, in accordance with the local policies where this device shall be used.

**Mui Scientific is not liable for any damage to the catheter, or harm to patients or personnel, caused by improper use of a disinfectant or procedure.**

REF Prefix	Description	REF Prefix	Description
CE1	Non-sleeve esophageal	LOSS	Lower Esophageal w/ Side holes & Sleeve
CE2	w/ Sleeve esophageal	UOSS	Upper Esophageal w/ Side holes & Sleeve
CNE2	w/ Sleeve Neonate Esophageal	PYS	Pyloric w/ sleeve
CPE1	Non-sleeve pediatric esophageal	ASS	Anorectal w/side holes & sleeve
CPE2	w/ Sleeve pediatric esophageal	APDSH	Antropyloroduodenal w/ side holes
CPR1	Non-sleeve pediatric anorectal	SISH	Small Intestines w/ side holes
CPR2	w/ Sleeve pediatric anorectal	OSH	Esophageal w/ side holes
CR1	Non-sleeve anorectal	ASH	Anorectal w/ side holes
CR2	w/ Sleeve anorectal		
CE4	Customized Esophageal	CE5	Customized Esophageal w/ sleeve
CR4	Customized Anorectal	CR5	Customized Anorectal w/ sleeve

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