| |
|---|
| DENTSLEEVE MANOMETRIC CATHETERS |
| INSTRUCTIONS FOR USE |
| Measure to be Taken Prior to Intubation |
| <ul style="list-style-type: none">Check that the design of the manometric catheter is specified for the intended measurements. Catheters are designed for recording from specific regions. Only use catheters that have been cleaned and sterilized according to instructions below. Inspect the catheter for any evidence of mechanical damage or imperfections. Do not use if damaged. Water fill channels that will not be used for measurements and seal catheter luer connectors with the luer plugs provided. |

Other equipment and Expertise Necessary for Measurements

- This catheter must be perfused with a manometric perfusion pump, at the recommended flow rate for achievement of technically valid measurements.
- Pressure amplification and recording equipment are also essential.
- This catheter must only be used under medical supervision by personnel who have been adequately trained in safe gastro-intestinal intubation, and in the performance of perfused manometric measurements.

Passage and Positioning of Catheter

- Dentsleeve upper gastro-intestinal manometric catheters are designed primarily for trans-nasal passage, after application of local anaesthetic.
- In the case of pharyngeal and oesophageal manometric catheters, correct positioning of the catheter is determined by evaluation of the pressure tracings.
- Special techniques are needed for reliable passage of manometric catheters across the pylorus and along the small intestine; verification of catheter position also requires special measures.

Positioning of Sleeve Catheters (if applicable)

- The Sleeve should be positioned so that it remains within the sphincter during its normal movement.
- After deeper insertion of the Sleeve, pull the catheter back into the sphincter, observing tracings from the Sleeve and side holes at either end of the Sleeve.
- Recording of a sphincteric pattern of pressure from a side hole at either end of the Sleeve indicates incorrect placement of the catheter relative to the sphincter.

Cleaning and Sterilisation

- This catheter is a re-usable device which has been validated as being able to withstand 50 cycles of steam autoclaving.
- If catheter has a balloon attached, plug the end of the luer connector to ensure no water enters into the balloon.
 - If catheter has impedance rings, do not immerse the Redel Connector in any liquid.
 - Immediately after extubation, immerse the catheter in a bowl of warm, mild detergent solution

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and wipe externally several times. Do not allow the uncleaned catheter to dry out.

- Flush the detergent solution gently through each channel with a 20ml syringe – Do not over-pressurise.
- Rinse the outside of the catheter in a bowl of clean water.
- Place the catheter on a clean towel, and cover it with a fold of the towel.
- Flush each channel first with water, then air – do not over-pressurise.
- Remove the luer connector plug from the balloon channel and autoclave with the catheter.
- Package the catheter for autoclaving, and autoclave without delay at 134°C for five minutes at 206 kPa or 30 psi, with a total cycle time of 30 minutes to allow for warm up and cool down time.

Intended Use

- Dentsleeve Manometric Catheters have been designed and manufactured for measurement of intralumal gastro-intestinal pressures in both adults and children when used with a manometric perfusion pump.
- This manometric catheter is to be used only under the supervision of a physician who has received adequate training in gastro-intestinal manometry.

Indications for Use

- Use of Dentsleeve Manometric Catheters is indicated when measurements of gastro-intestinal tract pressures are judged to be useful for determining management of patients with proven or suspected gastro-intestinal motor disorders.

Contraindications for Use

- Where there is a significantly increased risk of pulmonary aspiration associated with intubation which can not be adequately reduced by practical measures such as a longer than usual period of fasting prior to intubation.
- The presence of any stenosis, structural deformity or disease that results in an unacceptably increased risk of aspiration, perforation or any other adverse consequence arising from use of this assembly.

WARNING

- This devise is supplied non-sterile.
- Incorrect manometric perfusion pump operation can cause delivery of large volumes of gas through a manometric catheter during measurements.
- The use of "off-the-shelf" components from other manufacturers is not recommended. Only Dentsleeve components have been tested and validated to work safely and accurately with this catheter.
- Safe use of this device requires supervision by a physician trained in gastro-intestinal manometry.

USA only

Caution: US Federal law restricts this device to sale by or on the order of a physician.

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