## REUSABLE PVC MANOMETRIC CATHETERS INSTRUCTIONS FOR USE

### INTENDED PURPOSE

Reusable PVC Manometric Catheters are to be used in conjunction with manometric pumps and a computerized data processing system, to measure muscular contractions along the gastrointestinal system.

This manometric catheter is designed to be reusable, and to be used only under the supervision of a physician who has received adequate training in gastro-intestinal manometry.

### INDICATIONS FOR USE

Use of the Reusable PVC Manometric Catheters is indicated when measurements of gastro-intestinal tract pressure are judged to be useful for determining management of patients with proven or suspected gastro-intestinal motor disorders.

### CONTRAINDICATIONS

Where there is a significantly increased risk of pulmonary aspiration associated with intubation which cannot be adequately reduced by practical measures such as a longer than usual period of fasting prior to intubation.

The presence of any stenosis, structural deformity or disease that results in an unacceptably increased risk of aspiration, perforation or any other adverse consequence arising from use of the device.

**WARNINGS**  
Read through entire IFU to reduce any possible risks from misuse.

This device is supplied sterile.

Inspect catheter and packaging for any evidence of mechanical damage or imperfections. **Do not use if damaged.**

The catheter contains low amounts of phthalates. Use with caution for pregnant and pediatric cases. Reusable PVC Manometric Catheters must not be altered or modified in any way. Mui Scientific is not liable for personal injury and damage to property if original Mui Scientific parts are not being used.

Proper aseptic technique and universal barrier precautions (UBP) must apply.

Post-procedure manometric catheters are contaminated with body fluid. Immediately clean and disinfect according to the procedure stipulated in the accompanying Instructions For Use.

This catheter must only be used under medical supervision by personnel who have been adequately trained in safe gastro-intestinal intubation, and in the performance of perfused manometric measurements.

The detailed instructions contained within this IFU is sufficient for guidance of the usage of this equipment; no further training from the manufacturer is required.

Do not use devices beyond their prescribed lifetime or validated number of uses; injuries may result. Catheter is to be used in the anatomy and for the age group for which is intended, as dictated by the device name/description on the label.

### TARGET POPULATION

There are no specific intended patient populations.

The use for specific patients is at the discretion of the medical personnel who have received professional training in gastrointestinal manometry.

### CLINICAL BENEFIT

Manometric Catheters are an integral part of a motility system, for the measurement and assessment of the motor function of the gastrointestinal tract in both, adults and children. The clinical benefit of the use of the Reusable PVC Manometric Catheters, and this overall technique, is to enable the characterization of gastrointestinal patterns in the healthy and diseased patients, thus allowing the detection of motility disorders and the assessment of the therapy for patients suffering with problems along the gastrointestinal tract.

### INSTRUCTIONS FOR USE

Catheters are to be connected to pressure transducers, which are to be connected to a manometric pump and a computerized data processing system. Together as a system, these devices are used to measure muscular contractions along the gastrointestinal system. The type of pressure transducers and manometric pump used are dictated by the computerized data processing system being used for the motility study.

The Reusable PVC Manometric Catheters are compatible with any pressure transducer that has a standard Luer lock connector. Labeled female luer connector ends of the catheter are to be connected to the corresponding pressure transducer's male Luer. A male/female luer connection adaptor may be used, if required.

Catheters of the region of the gastrointestinal tract to be assessed corresponds with the variant type of the catheter device (esophageal, gastric, and small bowel catheters are to be intubated through the nasal or oral cavity; anorectal and colonic catheters are to be intubated through the rectum).

Water-based lubricants may be used to assist in insertion of catheter. **Do not use silicone-based lubricants on balloons.**

For anorectal catheters with a balloon, the balloon can be filled with air, to assess the patient's rectal sensations. The balloon may also be expelled by the patient, to assess the neuromuscular function of the rectum.

Post-procedure catheters are contaminated with body fluid; only re-use catheters that have been cleaned and reprocessed according to instructions below.

Do not use silicone-based lubricants on balloons.

The recommended perfusion flow rate for catheters with 1-12 channels is 0.6mL/min. Do not exceed flow rates greater than 6mL/min.

The recommended perfusion flow rate for catheters with 13+ channels is 0.15mL/min. Do not exceed flow rates greater than 1.5mL/min.

Polyisoprene (compliant) balloons – inflate up to rated volume; do not exceed double the rated volume; risk of balloon bursting.

Barostat (non-compliant) balloons – inflate up to rated volume; do not exceed rated volume; do not exceed pressures of 60mmHg; risk of balloon bursting.

### CHEMICAL CLEANING / DISINFECTING INSTRUCTIONS

The following reprocessing procedures have been validated by the manufacturer for the user to follow within their own facilities.

#### 1. Cleaning

- If catheter has a balloon attached, plug the end of the luer connector to ensure no liquid enters into the balloon.
- Wash catheter in warm water with enzymatic detergent.
- Flush detergent through all channels including unused channels. Remove all foreign matter. Wipe and clean external surface.
- Rinse catheter and flush all channels with clean fresh water to remove detergent. Repeat as needed.
- Blow-dry all channels with air and wipe excessive water from external surface.

#### 2. Chemical Disinfection (high-level disinfection)

- Flush disinfectant through all channels including unused channels. Ensure the entire catheter is submerged in disinfectant and all channels are filled with disinfectant.
- Soak catheter in disinfectant solution according to solution manufacturer's instructions.

#### Caution:

- The catheter contains Polyvinyl Chloride (PVC), Polycarbonate, Stainless Steel, Urethane Acrylate oligomer, N-Dimethylacrylamide and Isobornyl Acrylate. If a balloon is attached it also contains Polyisoprene and Silk Suture.
- Check with disinfectant manufacturer to ensure their disinfectant can be used safely on the materials listed above.
- Sporox (Sultan Chemicals Inc.), Korsolex Extra (Bode Chemie GmbH & Co.), OPTIM CS (SciCan Ltd.) and Adaspor Plus Ready to Use (IMS Group) are compatible with the catheter and should be used according to the solution's manufacturer's instructions.
- DO NOT USE CIDEX PRODUCTS TO DISINFECT ANY CATHETER THAT HAS A BALLOON.**
- Do not soak the catheter in alcohol.
- Do not steam autoclave the catheter.

#### 3. Final Rinsing and Blow Drying

- Rinse catheter and flush each channel with sterile water to remove disinfectant. Repeat as needed.
- Blow-dry all channels with air and wipe all water from external surface.
- The following reprocessing procedures have been validated by the manufacturer for the user to follow within their own facilities.

#### STORAGE

Ensure all channels and external surfaces are completely dry; store catheter dry in open air. **Do not store in a sealed container.**

Temperature limit:  
15°C - 30°C (60.8°F - 86.0°F)

Humidity limit:  
10% - 75%

Lifetime of product – 5 years or 50 uses whichever occurs first.

Time period for each use – 30 min inside the gastrointestinal system.

It is the responsibility of the user to manually track the age and number of uses of the device (User may contact Mui Scientific for a sample tracking form).

**Manufacturer's Warranty** – 6 months from Manufacture Date for any manufacturer's defects.

For each catheter system should report any serious incident that has occurred in relation to this device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Dispose of the device safely, in accordance with the local policies where this device shall be used. **Mui Scientific is not liable for any damage to the catheter, or harm to patients or personnel, caused by improper use of a disinfectant or procedure.**

Electronic copy of Instructions For Use is available on Mui Scientific website: [www.muiscientific.com](http://www.muiscientific.com)

| REF Prefix | Description          | REF Prefix | Description                    |
|------------|----------------------|------------|--------------------------------|
| R          | Esophageal           | C7         | Customized Reusable Esophageal |
| S          | Small Bowel          | CB-CG      | Barostat Gastric               |
| PE         | Pediatric Esophageal | CB-CR      | Barostat Anorectal             |
| PR         | Pediatric Anorectal  |            |                                |

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