REUSABLE PVC MANOMETRIC CATHETERS INSTRUCTIONS FOR USE

INTENDED PURPOSE

Reusable PVC Manometric Catheters are to be used in conjunction with manometric pumps and a computerized data processing system, to measure muscular contractions along the gastrointestinal system.

This manometric catheter is designed to be reusable, and to be used only under the supervision of a physician who has received adequate training in gastro-intestinal manometry.

INDICATIONS FOR USE

Use of the Reusable PVC Manometric Catheters is indicated when measurements of gastro-intestinal tract pressure are judged to be useful for determining management of patients with proven or suspected gastro-intestinal motor disorders.

CONTRAINDICATIONS

Where there is a significantly increased risk of pulmonary aspiration associated with intubation which cannot be adequately reduced by practical measures such as a longer than usual period of fasting prior to intubation.

The presence of any stenosis, structural deformity or disease that results in an unacceptably increased risk of aspiration, perforation or any other adverse consequence arising from use of the assembly.

WARNINGS

This device is supplied non-sterile.

Proper aseptic technique and universal barrier precautions (UBP) must apply.

Post-procedure manometric catheters are contaminated with body fluid.

Immediately clean and disinfect according to the procedure stipulated in the accompanying Instructions For Use.

This catheter must only be used under medical supervision by personnel who have been adequately trained in safe gastro-intestinal intubation, and in the performance of perfused manometric measurements.

INSTRUCTIONS FOR USE

- Inspect catheter prior to use; do not use if damaged.
- Catheters are to be used with a manometric pump and a computerized data processing system, to measure muscular contractions along the gastrointestinal system.
- Labeled luer connector ends are to be connected to the corresponding channel on the manometric pump.
- Post-procedure catheters are contaminated with body fluid; only re-use catheters that have been cleaned and reprocessed according to instructions below.
- Do not use silicone-based lubricants on balloons.

CHEMICAL CLEANING / DISINFECTING INSTRUCTIONS

1. Cleaning

- If catheter has a balloon attached, plug the end of the luer connector to ensure no liquid enters into the balloon.
- Wash catheter in warm water with enzymatic detergent.
- Flush detergent through all channels including unused channels. Remove all foreign matter. Wipe and clean external surface.
- Rinse catheter and flush all channels with clean fresh water to remove detergent. Repeat as needed.
- Blow-dry all channels with air and wipe excessive water from external surface.

2. Chemical Disinfection (high-level disinfection)

- Flush disinfectant through all channels including unused channels. Ensure the entire catheter is submerged in disinfectant and all channels are filled with disinfectant.
- Soak catheter in disinfectant solution according to solution manufacturer's instructions.

Caution:

- The catheter contains Polyvinyl Chloride (PVC), Polycarbonate, Stainless Steel (if a balloon is attached it also contains Polyisoprene or Polyolefin, and Silk Sutures). Check with disinfectant manufacturer to ensure their disinfectant can be used safely on the materials listed above.
- Cidex-OPA (Johnson & Johnson), Sporex (Sultan Chemists Inc.), Korsolex Extra (Bode Chemie GmbH & Co.), OPTIM CS (SciCan Ltd.) and Adaspor Plus Ready to Use (IMS Group) are compatible with the catheter and should be used according to the solution's manufacturer's instructions.
- DO NOT USE CIDEX PRODUCTS TO DISINFECT ANY CATHETER THAT HAS A BALLOON.
- Do not soak the catheter in alcohol.
- Do not steam autoclave the catheter.

3. Final Rinsing and Blow Drying

- Rinse catheter and flush each channel with sterile water to remove disinfectant. Repeat as needed.
- Blow-dry all channels with air and wipe all water from external surface.

ETO GAS REPROCESSING

- Ensure pigtail luer ends are open (remove all luer caps) prior to reprocessing.
- Reprocess at a temperature of 55°C for a minimum full-cycle of 180 minutes.
- Aerate at a temperature of 51-59°C for a minimum of 12 hours.

STORAGE

Ensure all channels and external surfaces are completely dry; store catheter dry in open air. Do not store in a sealed container.

Temperature limit:
15°C - 30°C (60.8°F - 86.6°F)

Humidity limit:
10% - 75%

Lifetime of product – 5 years or 50 uses whichever occurs first.

Manufacturer's Warranty – 6 months from Manufactured Date for any manufacturer's defects.

User is to contact the manufacturer or competent authority to report serious incidents involving this device.

Dispose of the device safely, in accordance with the local policies where this device shall be used.

Mui Scientific is not liable for any damage to the catheter, or harm to patients or personnel, caused by improper use of a disinfectant or procedure.

REF Prefix	Description	REF Prefix	Description
E	Esophageal	PR	Pediatric Rectal
R	Rectal	C7	Customized Reusable
S	Small Bowel	CB-CE	Barostat Esophageal
PH7	Acid Loading	CB-CG	Barostat Gastric
PE	Pediatric Esophageal	CB-CR	Barostat Rectal

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