## SINGLE-USE PVC MANOMETRIC CATHETER (ANORECTAL BALLOON EXPULSION – AIR ONLY) INSTRUCTIONS FOR USE

### INTENDED PURPOSE

The Anorectal Balloon Expulsion Catheter is a specific Single-Use PVC Manometric Catheter model which isolates the anorectal balloon function of an anorectal manometric catheter. The Anorectal Balloon Expulsion Catheter is used to evaluate the muscular contractions of the rectum.

### INDICATIONS FOR USE

Use of the Anorectal Balloon Expulsion Catheter is indicated when there is a need or suspicion of an anorectal disorder.

### CONTRAINDICATIONS

The presence of any stenosis, structural deformity or disease within the rectum, that results in an unacceptably increased risk of distention, perforation, or any other adverse consequence arising from use of this device.

### WARNINGS

Read through entire IFU to reduce any possible risks from misuse.

This catheter is supplied non-sterile.

Inspect catheter and packaging for any evidence of mechanical damage or imperfections. **Do not use if damaged.**

This catheter must not be altered or modified in any way. Mui Scientific is not liable for personal injury and damage to property if original Mui Scientific parts are not being used.

Proper aseptic technique and universal barrier precautions (UBP) must apply.

**This catheter is designed for single-use only; do not insert into patient more than once.**

Attempts to reprocess, re-sterilize, and/or re-use may lead to failure of the device and/or transmission of disease.

This catheter must only be used under medical supervision by personnel who have been adequately trained.

The detailed instructions contained within this IFU is sufficient for guidance of the use of this equipment; no further training from the manufacturer is required.

Do not use devices beyond their prescribed lifetime or shelf life; injuries may result. This catheter is to be used in the anatomy and for the age group for which is intended, as dictated by the device name/description on the label.

**Inflate balloon with air only;** risk of balloon detachment or rupture if balloon filled with water. (Use catheter SR2B if you wish to inflate balloon with water).

### TARGET POPULATION

There are no specific intended patient populations.

The use for specific patients is at the discretion of the medical personnel who have received professional training in gastrointestinal manometry.

### CLINICAL BENEFIT

The Anorectal Balloon Expulsion Catheter is to be inserted anorectally, for pressure measurements of the rectum. The balloon is to be filled with air to stimulate the rectum.

### BALLOON EXPULSION FUNCTION- INSTRUCTIONS FOR USE

1. Water-based lubricants may be used to assist in insertion of catheter. **Do not use silicone-based lubricants on balloons.**
2. Insert the catheter through the anus to a depth of 8-10cm.
3. Attach a 3-way stopcock to the catheter and fill the balloon with 30-60ml of air (do not fill with water; balloon detachment would occur).
4. Close off the balloon port, and remove syringe.
5. Have the patient go to the toilet / commode and attempt to expel the catheter.
6. Expulsion of the catheter within 2-5 minutes indicates normal function.
7. Post-procedure catheters are contaminated with body fluid, and must be immediately dispose of the catheter to avoid re-use.

### BIO-FEEDBACK FUNCTION – INSTRUCTIONS FOR USE

For handheld bio-feedback devices:

1. With a 3-way stopcock, attach the catheter to the male luer end of the stopcock, a syringe to one female luer end and an extension tube (with a male and female luer end) to the other female luer end of the stopcock.
2. Water-based lubricants may be used to assist in insertion of catheter. **Do not use silicone-based lubricants on balloons.**
3. Insert the catheter through the anus to a depth of 8-10cm.
4. Attach the other end of the extension tube to the handheld device and fill the balloon with air accordingly (30-60ml). Close off the syringe port, and remove syringe.
5. Any squeeze pressure by the patient will be displayed on the handheld device.
6. Post-procedure catheters are contaminated with body fluid, and must be immediately dispose of the catheter to avoid re-use.

### STORAGE

Store catheter in a dry and clean location.

Temperature limit:  
15°C - 30°C (60.8°F - 86.6°F)

Humidity limit:  
10% - 75%

### Shelflife of product – 2 years.

**Lifetime – 30 min inside the gastrointestinal system of the human body.**

**Manufacturer's Warranty** – 6 months from date of purchase from Mui Scientific for any manufacturer's defects.

The user and/or patient should report any serious incident that has occurred in relation to this device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Dispose of the device safely, in accordance with the local policies where this device shall be used.

**Mui Scientific is not liable for any damage to the catheter, or harm to patients or personnel, caused by improper use.**

Electronic copy of Instructions For Use is available on Mui Scientific's website: [www.muiscientific.com](http://www.muiscientific.com)

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Authorizes Representative in EU: Advena Ltd.
Authorizes Representative in Switzerland: Arazy Group Swiss GmbH
Authorizes Representative in UK: Advena Ltd.