# Mui Scientific

Division of  
H&A Mui  
Enterprises Inc.

145 Traders Blvd. E., Unit #34, Mississauga, ON, Canada, L4Z 3L3  
Tel: (905) 890-5525 Toll Free: 1-800-303-6611 Fax: (905) 890-3523  
Website: [www.muiscientific.com](http://www.muiscientific.com) Email: [mail@muiscientific.com](mailto:mail@muiscientific.com)

## SINGLE-USE PVC MANOMETRIC CATHETERS

### INSTRUCTIONS FOR USE

**Avoid inserting/using Single Use catheter into any patient more than once INTENDED PURPOSE**

Single-Use PVC Manometric Catheters are to be used in conjunction with manometric pumps and a computerized data processing system, to measure muscular contractions along the gastrointestinal system. This manometric catheter is designed for single-use, and is to be used only under the supervision of a physician who has received adequate training in gastro-intestinal manometry.

### INDICATIONS FOR USE

Use of the Single-Use PVC Manometric Catheters is indicated when measurements of gastrointestinal tract pressures are judged to be useful for determining management of patients with proven or suspected gastro-intestinal motor disorder.

### CONTRAINDICATIONS

Where there is a significantly increased risk of pulmonary aspiration associated with intubation which cannot be adequately reduced by practical measures such as a longer than usual period of fasting prior intubation. The presence of any stenosis, structural deformity or disease that results in an unacceptably increased risk of aspiration, perforation or any other adverse consequence arising from use of this assembly.

### WARNINGS

Read through entire IFU to reduce any possible risks from misuse.

This device is supplied non-sterile.

The catheter contains low amounts of phthalates. Use with caution for pregnant and pediatric cases.

Single-Use PVC Manometric Catheters must not be altered or modified in any way. Mui Scientific is not liable for personal injury and damage to property if original Mui Scientific parts are not being used.

Proper aseptic technique and universal barrier precautions (UBP) must apply. The Single-Use Manometric Catheters are designed for single-use only. Attempts to reprocess, re-sterilize, and/or re-use may lead to failure of the device and/or cross-contamination.

Post-procedure manometric catheters are contaminated with body fluid.

Immediately dispose of device to prevent re-use.

This catheter must only be used under medical supervision by personnel who have been adequately trained in safe gastro-intestinal intubation, and in the performance of perfused manometric measurements.

The detailed instructions contained within this IFU is sufficient for guidance of the usage of this equipment; no further training from the manufacturer is required. Do not use devices beyond their prescribed lifetime.

Catheter is to be used in the anatomy for which it is intended, as dictated by the device name/description on the label.

Catheter is to be used for the age group for which it is intended, as dictated by the device name/description on the label.

### TARGET POPULATION

There are no specific intended patient populations.

The use for specific patients is at the discretion of the medical personnel who have received professional training in gastrointestinal manometry.

### CLINICAL BENEFIT

Manometric Catheters are an integral part of motility system, for the measurement and assessment of the motor function of the gastrointestinal tract in both, adults and children. The clinical benefit of the Single Use PVC Manometric Catheters, and this overall technique, is to enable the characterization of gastrointestinal motor patterns in healthy and diseased patients, thus allowing the detection of motility disorders and the assessment of therapy suffering with problems along their gastrointestinal tract.

### INSTRUCTIONS FOR USE

Inspect catheter and packaging for any evidence of mechanical damage or imperfections.

### DO NOT USE IF DAMAGED.

Catheters are to be connected to pressure transducers, which are to be connected to a manometric pump and a computerized data processing system. Together as a system, these devices are used to measure muscular contractions along the gastrointestinal system. The type of pressure transducers and manometric pump used are dictated by the computerized data processing system being used for the motility study.

The Single Use PVC Manometric Catheters are compatible with any pressure transducer that has a standard luer lock connector. Labeled female luer connector ends of the catheter are to be connected to the corresponding pressure transducer's male luer. A male/female luer connection adaptor may be used, if required.

Ensure the region of the gastrointestinal tract to be assessed corresponds with the variant type of the catheter device.

Esophageal, gastric, and small bowel catheters are to be intubated through the nasal or oral cavity.

Anorectal catheters are to be intubated through the rectum.

For anorectal catheters with a balloon, the balloon can be filled with water or air, to assess the patient's rectal sensations. The balloon may also be expelled by the patient, to assess the neuromuscular function of the rectum.

Post-procedure manometric catheters are contaminated with body fluid. Immediately dispose of device to prevent re-use.

Do not use silicone-based lubricants on balloons.

The recommended perfusion flow rate for catheters with 1-12 channels is 0.6mL/min.

The recommended perfusion flow rate for catheters with 13+ channels is 0.15mL/min.

Polyisoprene / Polyolefin (compliant) balloons – inflate up to rated volume; do not exceed double the rated volume; risk of balloon bursting.

### STORAGE

Store catheter in a dry and clean location.

Temperature limit:  
15°C - 30°C (60.8°F - 86.6°F)

Humidity limit:  
10% - 75%

**Shelf life: 2 years for catheter with balloon and 3 years for catheter without balloon.**

**Lifetime: 30 min inside the gastrointestinal system of the human body.**

The user and/or patient should report any serious incident that has occurred in relation to this device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Dispose of the device safely, in accordance with the local policies where this device shall be used.

**Mui Scientific is not liable for any damage to the catheter, or harm to patients or personnel, caused by improper use.**

Electronic copy of Instructions For Use is available on Mui Scientific website:

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