## SINGLE-USE PVC MANOMETRIC CATHETERS INSTRUCTIONS FOR USE

### INTENDED PURPOSE

Single-Use PVC Manometric Catheters are to be used in conjunction with manometric pumps and a computerized data processing system, to measure muscular contractions along the gastrointestinal system. This manometric catheter is designed for single-use, and is to be used only under the supervision of a physician who has received adequate training in gastro-intestinal manometry.

### INDICATIONS FOR USE

Use of the Single-Use PVC Manometric Catheters is indicated when measurements of gastrointestinal tract pressures are judged to be useful for determining management of patients with proven or suspected gastro-intestinal motor disorder.

### CONTRAINDICATIONS

Where there is a significantly increased risk of pulmonary aspiration associated with intubation which cannot be adequately reduced by practical measures such as a longer than usual period of fasting prior intubation.

The presence of any stenosis, structural deformity or disease that results in an unacceptably increased risk of aspiration, perforation or any other adverse consequence arising from use of this assembly.

### WARNINGS

This device is supplied non-sterile.

Proper aseptic technique and universal barrier precautions (UBP) must apply.

The Single-Use Manometric Catheters are designed for single-use only. Attempts to reprocess, restirile, and/ or re-use may lead to failure of the device and/ or cross-contamination.

Post-procedure manometric catheters are contaminated with body fluid. Immediately dispose of device to prevent re-use.

This catheter must only be used under medical supervision by personnel who have been adequately trained in safe gastro-intestinal intubation, and in the performance of perfused manometric measurements.

### INSTRUCTIONS FOR USE

Inspect catheter and packaging for any evidence of mechanical damage or imperfections. **DO NOT USE IF DAMAGED.** Catheters are to be used with a manometric pump and a computerized data processing system, to measure muscular contractions along the gastrointestinal system.

Labeled luer connector ends are to be connected to the corresponding channel on the manometric pump.

Post-procedure manometric catheters are contaminated with body fluid. Immediately dispose of device to prevent re-use.

Do not use silicone-based lubricants on balloons.

### STORAGE

Store catheter in a dry and clean location.

Temperature limit:  
15°C - 30°C (60.8°F - 80.6°F)

Humidity limit:  
10% - 75%

**Shelf life and lifetime of product – 2 years for catheter with balloon and 3 years for catheter without balloon.**

User is to contact the manufacturer or competent authority to report serious incidents involving this device.

Dispose of the device safely, in accordance with the local policies where this device shall be used.

**Mui Scientific is not liable for any damage to the catheter, or harm to patients or personnel, caused by improper use.**

REF Prefix	Description
SE	Single-Use Esophageal
SR	Single-Use Rectal
SPE	Single-Use Pediatric Esophageal
SPRB	Single-Use Pediatric Rectal w/ Balloon
S7	Customized Single-Use

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