EN: DENTSLEEVE SILICONE MANOMETRIC CATHETERS
INSTRUCTIONS FOR USE

INTENDED PURPOSE

Dentsleeve Silicone Manometric Catheters are to be used in conjunction with a manometric pump and a computerized data processing system, to measure muscular contractions along the gastrointestinal system.

This manometric catheter is designed to be reusable, and is to be used only under the supervision of a physician who has received adequate training in gastro-intestinal manometry.

INDICATIONS FOR USE

Use of this catheter is indicated when measurements of gastro-intestinal tract pressures are judged to be useful for determining management of patients with proven or suspected gastro-intestinal motor disorders.

CONTRAINDICATIONS

Where there is a significantly increased risk of pulmonary aspiration associated with intubation which cannot be adequately reduced by practical measures such as a longer than usual period of fasting prior to intubation.

The presence of any stenosis, structural deformity or disease that results in an unacceptable increased risk of aspiration, perforation or any other adverse consequence arising use of this device.

WARNINGS

Read through entire IFU to reduce any possible risks from misuse.

This device is supplied non-sterile.

Inspect catheter and packaging for any evidence of mechanical damage or imperfections. **Do not use if damaged.**

Dentsleeve Manometric Catheters must not be altered or modified in any way. Mui Scientific is not liable for personal injury and damage to property if original Mui Scientific parts are not being used. Proper aseptic technique and universal precautions (UBP) must apply.

Post-procedure manometric catheters are contaminated with body fluid. Immediately clean and disinfect according to the procedures stipulated in the accompanying Instructions for Use.

This catheter must only be used under medical supervision by personnel who have been adequately trained in safe gastro-intestinal intubation, and in the performance of perfused manometric measurements.

The detailed instructions contained within this IFU is sufficient for guidance of the usage of this equipment, no further training from the manufacturer is required.

Do not use devices beyond their prescribed lifetime or validated number of uses; injuries may result. Use of this device is not to be used in the anatomy for the age group for which it is intended, as dictated by the device name/description on the label.

TARGET POPULATION

There are no specific intended patient populations.

The use for specific patients is at the discretion of the medical personnel who have received professional training in gastrointestinal manometry.

CLINICAL BENEFIT

Manometric Catheters are an integral part of a motility system, for the measurement and assessment of the motor function of the gastrointestinal tract in both, adults and children. The clinical benefit of the Dentsleeve Silicone Manometric Catheters, and the overall technique, is to enable the characterization of gastrointestinal motor patterns in healthy and diseased patients, thus allowing the detection of motility disorders and the assessment of therapy for patients with problems along their gastrointestinal tract.

INSTRUCTIONS FOR USE


- Water fill channels that will not be used for measurements and seal catheter luer connector with the luer plugs provided.
 - Labeled luer connector ends are to be connected to the corresponding channel on the manometric pump.
 - Ensure the region of the gastrointestinal tract to be assessed corresponds with the variant type of the catheter device (esophageal, gastric, and small bowel catheters are to be intubated through the nasal or oral cavity; anorectal and colonic catheters are to be intubated through the rectum).
 - The recommended perfusion flow rate for catheters with 1-12 channels is 0.6mL/min. Do not exceed flow rates greater than 6mL/min.
 - The recommended perfusion flow rate for catheters with 13+ channels is 0.15mL/min. Do not exceed flow rates greater than 1.5mL/min.
 - For anorectal catheters with a balloon, the balloon can be filled with air only, to assess the patients rectal sensations. The balloon may also be expelled by the patient, to assess the neuromuscular function of the rectum.
 - Post-procedure catheters are contaminated with body fluid; only re-use catheters that have been cleaned and reprocessed according to instructions below.
 - Proper aseptic technique and universal barrier precautions (UBP) must apply.
- POSITIONING OF SLEEVE CATHETERS (IF APPLICABLE)**
- The Sleeve should be positioned so that it remains within the sphincter during its normal movement.
 - After deeper insertion of the Sleeve, pull the catheter back into the sphincter, observing tracings from the Sleeve and side holes at either end of the Sleeve.
 - Recording of a sphincter pattern of pressure from a side hole at either end of the Sleeve indicates placement of the catheter relative to the sphincter.

CLEANING AND REPROCESSING

- This catheter is a re-usable device which has been validated as being able to withstand 50 cycles of steam autoclaving.
- If catheter has balloon attached, plug the end of the luer connector to ensure no water enters into the balloon.
 - If catheter has impedance rings, do not immerse the Redel Connector in any liquid.
 - Immediately after expiration, immerse the catheter in a bowl of warm, mild enzymatic detergent and wipe externally several times. Do not allow the uncleaned catheter to dry out. Flush the detergent solution gently through each channel -- Do not over-pressurize.
 - Rinse the outside of the catheter in a bowl of clean water. Place the catheter on a clean towel, and cover it with a fold of the towel.
 - Flush each channel first with water, then air -- do not over-pressurize. Remove the luer connector plug from the balloon channel and autoclave with the catheter.
 - Package the catheter for autoclaving, and autoclave without delay at 134C for three minutes at 206kPa or 30psi with a total cycle time of 30 minutes to allow for warm up and cool down time. Dry catheter for 30 minutes.

STORAGE

- Ensure all channels and external surfaces are completely dry; store catheter dry in open air. Do not store in a sealed container.

Temperature limit:
15°C - 30°C (60.8°F - 86.6°F)

Humidity limit:
20% - 75%

Lifetime of product – 5 years or 50 uses whichever occurs first.

Time period for each use - 30 min inside the gastrointestinal system.

It is the responsibility of the user to manually track the age and number of uses of the device (User may contact Mui Scientific for a sample tracking form).

Manufacturer's Warranty – 6 months from Manufactured Date for any manufacturer's defects.

The user and/or patient should report any serious incident that has occurred in relation to this device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Dispose of the device safely, in accordance with the local policies where this device shall be used. **Mui Scientific is not liable for any damage to the catheter, or harm to patients or personnel, caused by improper use of a disinfectant or procedure.**

Electronic copy of Instructions For Use is available on Mui Scientific website: www.muiscientific.com

REF Prefix	Description	REF Prefix	Description
CE1	Non-sleeve esophageal	LOSS	Lower Esophageal w/ Side holes & Sleeve
CE2	w/ Sleeve esophageal	UOSS	Upper Esophageal w/ Side holes & Sleeve
CPE1	Non-sleeve pediatric esophageal	PYS	Pyloric w/ sleeve
CPE2	w/ Sleeve pediatric esophageal	ASS	Anorectal w/side holes & sleeve
CPR1	Non-sleeve pediatric anorectal	APDSH	Antropyloroduodenal w/ side holes
CPR2	w/ Sleeve pediatric anorectal	SISH	Small Intestines w/ side holes
CR1	Non-sleeve anorectal	OSH	Esophageal w/ side holes
CR2	w/ Sleeve anorectal	ASH	Anorectal w/ side holes
CE4	Customized Esophageal	CE5	Customized Esophageal w/ sleeve
CR4	Customized Anorectal	CR5	Customized Anorectal w/ sleeve

Basic UDI-DI: 0678467DENTSLEEVE9

Made in Canada

Authorizes Representative in UK: Advena Ltd.:

Pure Offices, Plato Close, Tachbrook Park, Warwick, CV34 6WE, UK

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