EN: PATCOM SINGLE-USE INTRODUCTOR INSTRUCTIONS FOR USE

INTENDED PURPOSE

The PatCom Single-Use Introducer is used for patient esophageal intubation to guide catheters and tubes via nasal or oral cavity. It is to be used in conjunction with endoscopes to allow visualization of the placement process.

INDICATIONS FOR USE

The PatCom Single-Use Introducer should be used for the transnasal or transoral placement of any catheter or tube into the esophagus that otherwise don't offer visualization of the placement process to the end user.

CONTRAINDICATIONS

When the outside diameter of the introducer is too large for the patients nasal or oral passages.

When the inside diameter of the introducer is too small for the outside diameter of the catheter/ tube to be placed.

Do not use in patients with esophageal bleeding and/ or lacerations; laryngeal perforations; trauma to teeth, gums and/ or pharynx, aspiration pneumonia. Contraindications also include those specific to any endoscopic procedure, use of an overtube, or any endoscopic procedure performed with an overtube.

WARNINGS

Read through entire IFU to reduce any possible risks from misuse.

Inspect the device and packaging for any mechanical damage or imperfections. Do not use introducer if damaged.

The PatCom Single-Use Introducer is for transient use (under 24 hrs), and will come in direct contact with patient's mucosal lining.

The PatCom Single-Use Introducer is provided non-sterile; proper aseptic technique and universal barrier precautions (UBP) must apply.

The device is intended to be non-reusable. Attempts to reprocess, re-sterilize and/or re-use may lead to failure of the device and/or transmission of disease. Post-procedure introducers are contaminated with body fluid; immediately disposed of device to prevent re-use.

The PatCom Single-Use Introducer is to be used in a hospital or clinical setting and under the supervision of a qualified healthcare professional who has received professional training in using the equipment.

Use of over tubes have been associated with perforation. Such complications should be mentioned during and after the procedure and treated appropriately, if detected. Residual water within the body of the introducer may pose an aspiration risk if not removed prior to intruder extubation.

Do not use devices beyond their prescribed lifetime.

Catheter is to be used in the anatomy for which it is intended, as dictated by the device name/description on the label.

Catheter is to be used for the age group for which it is intended, as dictated by the device name/description on the label.

TARGET POPULATION

There are no specific intended patient populations.

The use for specific patients is at the discretion of the medical personnel who have received professional training in the gastrointestinal system.

CLINICAL BENEFIT

Single-Use Introducer provides guidance of medical devices through the nose or mouth into the pharynx and esophagus. It assists in avoiding misplacement of catheters and probes into the undesired location like the trachea.

INSTRUCTIONS FOR USE

1. Prepare clean working area.
2. Choose an Introducer with an inner (ID) that is at least 0.4mm larger than the outer diameter (OD) of the endoscope and/or catheter that is to be placed inside the Introducer.
3. Place introducer, endoscope, catheter or tube, and lubricant on working area.
4. Apply lubricant generously to outside of Introducer (first 5 cm from tip), outside of endoscope (first 5 cm from tip), and outside of catheter or tube (first 5 cm from tip).
5. Slide Introducer starting with the flared end over insertion portion of endoscope as far up as possible.
6. For nasal intubation:
 - a) While holding Introducer and endoscope, insert both into the nasal cavity of the patient and slide forward utilizing the visualization capabilities of the endoscope.
 - b) When the post pharyngeal wall is reached, utilize the angulated tip of the endoscope to navigate introducer and endoscope into the pharynx.
 - c) While instructing the patient to swallow, move Introducer/ endoscope forward into the esophagus until the flared end of the introducer is at the tip of the patient's nose.
7. For oral intubation:
 - a) While holding the Introducer and endoscope, insert both into the oral cavity of the patient and slide forward utilizing the visualization capabilities of the endoscope.
 - b) Slide endoscope and introducer laterally over the epiglottis until the larynx is visible.
 - c) While instructing the patient to swallow, move the introducer and endoscope forward into the esophagus until the flared end of the introducer is at the opening of the patient's mouth.
8. While holding the Introducer at the flared end in place at the patient's nose/ mouth, slowly pull out the endoscope until it is completely extracted.
9. While holding the Introducer at the flared end in place at the patient's nose/ mouth, slowly insert the lubricated catheter/ tube into the introducer and push through introducer to the desired position in the esophagus.
10. Slowly pull back the Introducer while maintaining position of catheter/ tube until the Introducer is fully extracted.
11. Hold flared end of Introducer with both hands and use tear away seams to peel introducer off the catheter/ tube. Dispose of introducer and begin desired test or procedure.

STORAGE

Store Introducer in dry and clean location.

Temperature limit:
15°C - 30°C (60.8°F - 86.6°F)

Humidity limit:
10% - 75%

Shelf life of product - 2 years. Lifetime of product - 30 minutes.

Manufacturer's Warranty - 6 months from Manufactured Date for any manufacturer's defects.

The user and/or patient should report any serious incident that has occurred in relation to this device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Dispose of the device safely, in accordance with the local policies where this device shall be used.

Mui Scientific is not liable for any damage to the introducer/ endoscope, or harm to patients or personnel, caused by improper use.

Electronic copy of Instructions For Use is available on Mui Scientific website: www.muiscientific.com

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