REF : S7-BR-1011  
: S7-BR-1018

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### EN: RAPID BAROSTAT BAG CATHETER INSTRUCTIONS FOR USE

#### INTENDED PURPOSE

The Rapid Barostat Bag (RBB) Catheter is designed specifically to be used with the Rapid Barostat Bag Pump to measure rectal capacity as well as rectal volume at points of sensation, urge, and discomfort, to determine hypersensitivity, hyposensitivity, or normal rectal sensations. These determinations aid in the diagnosis of anorectal disorders, such as fecal incontinence and constipation.

#### INDICATIONS FOR USE

Use of the RBB Catheter (together with RBB Pump) is indicated when there is a need or suspicion of an anorectal disorder or for evaluation of other rectal motor sensory functions.

#### CONTRAINDICATIONS

The presence of any stenosis, structural deformity or disease within the rectum, that results in an appreciably increased risk of distention, perforation, or any other adverse consequence arising from use of this device.

#### WARNINGS

Read through entire IFU to reduce any possible risks from misuse.

This device is supplied non-sterile. Inspect catheter and packaging for any evidence of mechanical damage or imperfections. **Do not use if damaged.**

The catheter must not be altered or modified in any way. Mui Scientific is not liable for personal injury and damage to property if original Mui Scientific parts are not being used. Proper aseptic technique and universal barrier precautions (UBP) must apply.

**The catheter is designed for single use only; do not insert into patient more than once.** Attempts to reprocess, sterilize, and/or reuse may lead to failure of the device and/or cross-contamination. The RBB Catheter (together with RBB Pump) must be used in a medical setting and should be operated by a qualified medical professional trained in gastrointestinal manometry.

Digital rectal examination should be performed prior to inserting the catheter into the patient.

Do not use the device on patients with pre-existing rectal pain or bleeding (at the discretion of the physician).

Do not use catheter beyond their prescribed lifetime or shelf life; injuries may result.

This catheter is to be used in the anatomy and for the age group which is intended, as dictated by the device name/description on the label.

#### TARGET POPULATION

There are no specific intended patient populations. The use for specific patients is at the discretion of the medical personnel who have received professional training in gastrointestinal manometry.

#### CLINICAL BENEFIT

The RBB Catheter is used together with RBB Pump for the inflation of a Barostat Bag to measure Rectal Capacity, Pressure, Compliance and Sensitivity. These findings can aid in the diagnosis of anorectal disorders such as fecal incontinence and constipation, as well as aid in surgical pre-screening, and determining rectal sensitivity of Inflammatory Bowel Disease and Irritable Bowel Syndrome patients.

#### INSTRUCTIONS FOR USE

Ensure RBB-pump has been turned ON, and the patient's information has been entered.

Have the patient assume the left lateral decubitus position.

Attach the luer end of the RBB catheter to the luer at the top left hand side of the RBB pump.

Press the 'Empty' button on the touch screen to remove any existing air from the RBB catheter.

Lubricate the RBB catheter and insert into the patients rectum until the 5cm marking is at the anus. Tape down the catheter to the buttocks.

When the patient has become acclimated, press the "Start" button to commence the "Capacity Measurement" study.

Follow the RBB pump prompts to commence the "Sensation Measurement" study.

Follow the RBB pump prompts to "Repeat" or "End of study". Upon completion of study, immediately dispose of catheter to prevent reuse and cross contamination.

#### STORAGE

Store catheter in a dry and clean location.

Temperature limit:  
15°C - 30°C (60.8°F - 86.6°F)

Humidity limit:  
10% - 75%

#### Shelf life of product – 2 years.

**Lifetime – 30 min inside the gastrointestinal system of the human body.**

**Manufacturer's Warranty** – 6 months from date of purchase from Mui Scientific for any manufacturer's defects.

**The user and/or patient should report any serious incident that has occurred in relation to this device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.**

Dispose of the device safely, in accordance with the local policies where this device shall be used.

**Mui Scientific is not liable for any damage to the catheter, or harm to patients or personnel, caused by improper use.** Electronic copy of Instructions For Use is available on Mui Scientific's website: [www.muiscientific.com](http://www.muiscientific.com)

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