





















PVC TUBING EXTRUSIONS

 Extrusion 8F9 8+1 Motility Lumina 0.5mm & 1.0mm 2.7mm	 Extrusion 10F4 4 Motility Lumina 0.8mm 3.2mm	 Extrusion 10F9 8+1 Motility Lumina 0.6mm & 1.25mm 3.2mm
 Extrusion 12F9 8+1 Motility Lumina 0.7mm & 1.3mm 3.9mm	 Extrusion 14F8 6+1+1 Motility Lumina 0.6mm, 1.3mm & 2.2mm 4.6mm	 Extrusion 14F9 8+1 Motility Lumina 0.8mm & 1.8mm 4.7mm
 Extrusion 14F13 12+1 Motility Lumina 0.6mm & 2.2mm 4.6mm	 Extrusion 12F25 24+1 Motility Lumina 0.3mm & 1.0mm 4.0mm	 Extrusion 14F25 24+1 Motility Lumina 0.43mm & 1.65mm 4.7mm

SILICONE TUBING EXTRUSIONS

 Extrusion E1 Standard 8+1 Lumina 0.74mm & 1.55mm 4.7mm	 Extrusion E4 Micro 4+1 Lumina 0.35 mm & 0.6mm 1.75mm	 Extrusion E6 Micro 9+1 Lumina 0.35mm & 0.75mm 2.0mm
 Extrusion E27 Standard 8+1 Lumina 0.65mm & 0.9mm 3.5mm	 Extrusion E30 Miniature 16+1 Lumina 0.5mm & 1.0mm 3.5mm	 Extrusion E55 Barostat Miniature 20+2+1 Lumina 0.4mm & 0.6mm Oval 1.9mm x 2.4mm 4.2mm
 Extrusion E56 Miniature 12+1 Lumina 0.4mm & 1.9mm 3.5mm	 Extrusion E57 Miniature 21+1 Lumina 0.4mm & 0.9mm 4.0mm	 Extrusion E63 Aspiration/Infusion Standard 2+1+1 Lumina 0.6mm & 0.8mm Oval 1.5mm x 1.2mm 3.0mm
 Extrusion E114 Standard 8+1 Lumina 0.67mm & 1.39mm 4.0mm	 Extrusion E157 Miniature 41 Lumina 0.4mm 4.0mm 5.2mm	