



America

# CERTIFICATE

No. U8 18 05 86815 004

Holder of Certificate: **Mui Scientific**  
**Division of H & A Mui Enterprises Inc.**

145 Traders Blvd. E, Unit 34  
Mississauga ON L4Z 3L3  
CANADA

Production Facility(ies): 86815

Certification Mark:



**C US**

Product: **General Medical Devices**  
**Manometric Pump**

Model(s): **Pumps: PIP-4-xxSS, P4-F-203, P4-F-210**  
**Accessories: PIP-5-2SS, PIP-5-3SS, PIP-6-xxSS**  
**Where xx is the number of channels and can be 4, 6, 8, 10, 12, 14, 16, 18, 20, 22, or 24.**

Parameters: Pump:  
Rated Input Voltage: 115 V~  
Rated Frequency: 60 Hz  
Rated Input Current: 2.0 A  
Duty Cycle: Max 5 minutes on, 115 minutes off  
Protection Class: I  
Degree of Protection: Type BF  
Accessories:  
Input Pressure: 117-373 kPa

Tested according to: CAN/CSA-C22.2 No. 60601-1:2014 excluding Biocompatibility (clause 11.7), Usability (clauses 12.2 and 15.1), and EMC (clause 17). ANSI/AAMI ES60601-1:2005/A1:2012-08 excluding Biocompatibility (clause 11.7), Usability (clauses 12.2 and 15.1), and EMC (clause 17).

The product was voluntarily tested according to the relevant safety requirements noted above. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited certification body.

Test report no.: 7169004280-000

Date, 2018-05-30

Page 1 of 1

