



America

CERTIFICATE

No. U8 15 11 86815 003

Holder of Certificate: **Mui Scientific
Division of H & A Mui Enterprises Inc.**

145 Traders Blvd. E
Mississauga ON L4Z 3L3
CANADA

Production Facility(ies): 86815

Certification Mark:



C US

Product: **General Medical Devices
Rapid Barostat Bag Pump**

Model(s): **P1-RBB-1**

Parameters:

Device:	
Rated Input Voltage:	5.0 Vdc
Rated Input Power:	2.5 A
Degree of Protection:	Type BF applied parts
Charger:	
Rated Input Voltage:	100-240 V~
Rated Frequency:	50-60 Hz
Rated Current:	1.0-0.5 A
Protection Class:	II

Tested according to: CAN/CSA C22.2 No. 60601-1:2014
excluding Biocompatibility (clause 11.7),
Usability (clauses 12.2, 15.1), and
EMC (clause 17).
ANSI/AAMI ES 60601-1:2005/R:2012
excluding Biocompatibility (clause 11.7),
Usability (clauses 12.2, 15.1), and
EMC (clause 17).

The product was voluntarily tested according to the relevant safety requirements noted above. It can be marked with the certification mark above. The mark must not be altered in anyway. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC Guide 67. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited certification body.

Test report no.: 240-1500009-000/KM1500009

Date, 2015-11-06

Page 1 of 2

