





Reusable Barostat Catheters

Type	Balloon Volume (mL)	Balloon Lumen ID (mm)	Measurement Lumen ID (mm)	Part Number
Esophageal	150	2.0	1.5	CB-CE-0001
Gastric	800	3.3	2.0	CB-CG-0001
Rectal	600	3.3	2.0	CB-CR-0001



Replacement Barostat Balloons (Pillow-Type)


(custom sizes available)



Type	Volume (mL)	Length (cm)	Diameter (cm)	Part Number
Esophageal	150	15	9	CT-BP150E
Anorectal	600	22	15	CT-BP600R
Gastric	800	25	16	CT-BP800G
Gastric	1000	25	18	CT-BP1000G

Replacement Barostat Balloons (Square-Type)

(custom sizes available)




Esophageal	150	13	6	CT-B150E
Anorectal	600	18	10	CT-B600R
Gastric	800	23	10	CT-B800G
Gastric	1000	25	10	CT-B1000G



Single-Use Barostat Catheters

Type	Balloon Volume (mL)	Balloon Lumen ID (mm)	Measurement Lumen ID (mm)	Part Number	
Esophageal	150	2.0	1.5	S7-CB-E001	(10 per box)
Gastric	800	3.3	2.0	S7-CB-G001	
Rectal	600	3.3	2.0	S7-CB-R002	



Single-Use Anorectal Balloon Expulsion / Biofeedback Catheter

	# of ch.	O.D. (mm)	Part Number	
	1	5.7	SR1B	(10 per box)

Single-Use Anorectal Rapid Barostat Bag (RBB) Catheter

	# of ch.	I.D. (mm)	Part Number	
	1	3.2	S7-BR-1011 (with female luer)	(10 per BOX)
	1	3.2	S7-BR-1018 (with male luer)	(10 per BOX)