



ANORECTAL BALLOON EXPULSION SERIES

Part Number	SR1B	SR2B	CR1B
Intended Use	The Anorectal Balloon Expulsion Catheter is used to evaluate the muscular contractions of the rectum.		
Reusable			 Up to 50 uses
Size of Balloon	Maximum Volume: 400cc		
Can be filled with		 or 	
Incremental Markings			
Cost	\$	\$\$	\$\$\$

For any inquiries or further details, please contact us: Tel: (905) 890-5525 | Email: mail@muiscientific.com